


Test result summary report

Results of sensitivity test performed on March 5, 2020



RNA concentration	Lot1	
	10^6 copies	Positive
10^5 copies	Positive	Positive
10^4 copies	Positive	Positive
10^3 copies	Positive	Positive
10^2 copies	Positive	Positive
Negative control	Negative	Negative

Samples of 10^6 copies, 10^5 copies, 10^4 copies, 10^3 copies, and 10^2 copies were prepared from the artificial RNA of the Nsp15 region of SARS-CoV2, and the detection sensitivity by the SmartAmp method was examined.

- 1) The sensitivity is to detect 10^2 copies (100 copies) in 30-minute measurement.
- 2) Negative control is negative even after 60 minutes.

The results are shown above.

March 6, 2020

Prepared by

RIKEN Preventive Medicine & Diagnosis Innovation Program (PMI) supervisor

何俊泰司

Test result manager

RIKEN Preventive Medicine & Diagnosis Innovation Program (PMI) Program director

Yoshihide Hayashizaki 林崎良英